


db lab supplies

AUTUMN OFFERS

valid until 22nd January 2010



FREE GIFTS

Order must be placed in one single transaction - one gift only per transaction. Gifts can not be used in conjunction with any other offer or discount scheme or gift campaign. Returns will be excluded from the order value. In the event of the order amount falling below the required gift order value you will be required to return the gift or pay for the free gift in full. Only applies to goods ordered from this leaflet!

spend over **£250**
and receive

Leowire
0.7mm 500g
worth **£39.72**



spend over **£500**
and receive

ix164 Adams
worth **£89.00**



spend over **1000**
and receive
ipod Nano
worth **£107.00**



NSK Ultimate 450 Bench Top Micromotor

Ultimate 450 is a laboratory micromotor with a bench-top control unit. Its microprocessor-controlled power offers versatility for every type of laboratory works. Providing high torque of 4.5 Ncm and wide speed range from 1,000 to 35,000 min⁻¹, Ultimate 450 is an ideal choice for all laboratory applications.

Features

- Optimal Control by Microprocessor
- Self-Diagnosis Function
- Smooth Operation
- Audible Check System
- Easy-to-operate Display Panel
- Maintenance Mode Setting (450/450K)
- Digital Speed Display (450/450K)
- Error Code Display Function (450/450K)

Ultimate 450 Systems



Contents:
Control Unit, Micromotor Handpiece,
Foot Control, Handpiece Stand

Specifications

Power Supply AC 230V 50/60Hz
Weight 2.6kg
Dimensions W172xD230xH94 (mm)



NRP £971.75
Autumn Offer
£870.00



PHONE **01535 650 570**



Orthodontic Solders

For all orthodontic soldering applications. Cadmium free. Silver content: 55%. Low melting point: 600°. Packs of 10g approx. LR0222-00 Wire form, 0.40mm in diameter LR0223-00 Tape form, 1 x 0.5mm LR0224-00 Rods, 1mm in diameter NRP £29.05

Special Offer
£24.70



Flux

In paste. For use with solders: LR0222-00, LR0223-00 and LR0224-00. Pack of 30g approx. LR0226-00

NRP £11.64

Special Offer
£9.90



Thermo Insulator Paste for Solders

It protects plaster and resin during soldering procedure. Excellent consistency. Easy to remove with water. Pack of 235ml. LR0227-01

NRP £24.01

Special Offer
£20.40

Stainless Steel Wire - Leowire


Chromium stainless steel alloy diamond drawn to obtain perfectly smooth, clean and round surface. Sold on a 500g reel.



- 0.5mm LC0450-05
- 0.6mm LC0450-06
- 0.7mm LC0450-07
- 0.8mm LC0450-08
- 0.9mm LC0450-09
- 1.0mm LC0450-10
- 1.1mm LC0450-11
- 1.2mm LC0450-12
- 1.3mm LC0450-13
- 1.4mm LC0450-14

£39.72/reel

Buy 3 GET 1 FREE



Quad Helix Arches

Made from 0.036"/0.9mm stainless steel with curved ends for easy insertion. No soldering needed. For maxillary and mandible with correct gingival clearance.

- X-small DB52-0000U
- Small DB52-0001U
- Medium DB52-0002U
- Large DB52-0003U

£6.00/each

Buy 4 GET 1 FREE



Ball Clasps

The perfect round shaped one piece ball clasp. Sold in packs of 100.

- 0.7mm DBL2-0007
- 0.8mm DBL2-0008
- 0.9mm DBL2-0009
- 1.0mm DBL2-0010

£14.30/pack

Buy 3 Packs GET 1 PACK FREE



Sheaths

Curved DB52-0010

NRP £30.60/pack 10

Special Offer
£27.50/pack 10



Adams Plier

IX164

The perfect plier for constructing perfect Adams Clasps. Max wire diameter 1.0mm

NRP £89.00

Special Offer
£79.00



Spring Former Plier

IX165

The perfect plier for the most accurate coils and springs. Max wire diameter 0.7mm

NRP £89.00

Special Offer
£79.00



Three Jaw Plier

IX701

The perfect plier for bending the arms of all rapid expansion screws. Max wire diameter 1.4mm

NRP £89.00

Special Offer
£75.65




Young Plier

IX703


The perfect plier of all labial arches requiring accurate "U" loops. Max wire diameter 0.8mm.

NRP £75.00

Special Offer
£75.00



The inner spring produces a constant force producing a controlled bio mechanical movement of the arch.



Spring Activated Expander

Ideal for the expansion of the upper arch. This spring expander incorporates a built-in stop to eliminate over-expansion. Reactivation is achieved with a special key usually once a month.

A LA2701-08

Light spring approx. 500g

B LA2702-08

Medium spring approx. 900g

£21.75/each

Buy 4 GET 1 FREE
on all expansion screws



← 7mm →
LA0805-11

← 9mm →
LA0805-14

Standard Expansion Screw

£16.36/pack of 10



← 7mm →
LA0803-11

← 9mm →
LA0803-14

Medium Expansion Screw

£16.36/pack of 10



← 6mm →
LA0800-10

← 8mm →
LA0800-12

Universal Expansion Screw

£17.85/pack of 10



Leocryl Orthodontic Self Curing Acrylic Resin

Designed especially for the construction of all orthodontics appliances. The special formula of polymer and monomer offers the following properties:

- Slow flow
- Accurate fit to plaster models
- Superior strength and hardness

The Leone Hydroflask or pressure pot is recommended for superior results eliminating porosity and residual monomer. Polymerisation is recommended at 40°C at a pressure of 2.5 Atm for 20 minutes.

Regular Kit

LR6222-95 - Clear LR6223-95 - Transparent Pink

The Kit includes; 1 Kg powder, 0.5lt liquid, spray flasks for powder/liquid, funnels, beakers and separating medium.

NRP £145.00

Special Offer
£123.00

Introductory Kit

LR6202-93 - Clear LR6203-93 - Transparent Pink

Kit includes; 60g powder, and 45ml liquid.

NRP £31.45

Special Offer
£25.00



Leone Hydroflask

Orthodontic cold curing acrylics require the polymerization under pressure technique. This kind of polymerization enables the vaporization of monomer and condensing of acrylic, which result in a perfect fit to the plaster model. Appliances are easily finished and can be immediately polished. The hydroflask allows you to make use of the polymerization under pressure technique easily and economically.


Technical Features

Inside dimensions:
100 x 115mm
Outside dimensions:
111 x 130mm
Pressure gauge
Exhaust Valve
Made of hard Aluminium
Built-in thermometer

LT1315-00

NRP £95.95

Special Offer
£81.55



Press for Hydroflask

Convenient press for Hydroflask T1310-00 and T1315-00, made of pressure cast aluminium. Extremely easy to use.

LT1311-00

NRP £79.95

Special Offer
£67.95



Acrylic Resin Coloured Liquids

Sold in 250ml bottles.

Red	LR6124-00
Purple	LR6125-00
Magenta	LR6126-00
Orange	LR6127-00
Yellow	LR6128-00
Green	LR6129-00
Turquoise	LR6130-00
Blue	LR6131-00

Also available in 0.5l bottles.

NRP £24.75

Special Offer
£19.95



Glitters

Attractive looking appliances are obtained by mixing a small quantity of glitters with acrylic powder. Glitters make appliances x-ray visible. Packs of 50g and a spray flask.

LR4300-00	Silver
LR4301-00	Gold
LR4302-00	Multicolour "Z"
LR4303-00	Blue
LR4304-00	Red

NRP £18.50

Special Offer
£15.00

db 
lab supplies

DB LAB SUPPLIES
UNIT 6 RYEFIELD WAY SILSDEN
WEST YORKSHIRE BD20 0EF
UNITED KINGDOM

TELEPHONE > +44[0] 1535 650 570
FAX > +44[0] 1535 656 969
EMAIL > sales@dblabssupplies.co.uk
WEB > www.dblabssupplies.co.uk

DB LAB SUPPLIES IS A DIVISION OF DB ORTHODONTICS LIMITED.