


Archshape

EuroLine™

Euroline is a European Arch Shape and the most popular in our range.
Available in Nickel Titanium, Thermal, Stainless Steel, Posted, Coated and Beta.




Triline™

Triline is a natural arch shape available in 3 sizes.


Features of Triline include:

- > No over expansion of the pre-molars.
- > Wider in the molar regions when compared to standard arch shapes.
- > Cuspid marking available on Stainless Steel.



Aesthetic Micro-Coated Archwires (Thermal NiTi, S.E. NiTi & Stainless Steel)

EuroLine™ European Arch Shape



- > Micro-Coated Aesthetic Archwires - only coated on labial surface.
 - > Micro-Coated Aesthetic Archwires - a huge advancement in coated wires provide an optimum cosmetic appearance.
 - > Micro-Coated Archwires will reduce friction and improve sliding mechanics compared to fully coated wire.
 - > Micro-Coated Archwires achieve a true archwire size as the PTFE coating is only .0005".
 - > The coating is resistant to staining and cracking even if the archwire shape is changed.
 - > All Micro-Coated Archwires are individually packaged & sealed in an autoclavable pack to protect your patients from cross infection.
 - > Micro-Coated labial surface provides an Optimum cosmetic Appearance.
 - > Pack size: 10
- ALL SIZES AVAILABLE

Nickel Titanium

Round	Upper	Lower
.012"	DB71-0012U	DB71-0012L
.014"	DB71-0014U	DB71-0014L
.016"	DB71-0016U	DB71-0016L
.018"	DB71-0018U	DB71-0018L
.020"	DB71-0020U	DB71-0020L

Stainless Steel

Round	Upper	Lower
.012"	DB72-0012U	DB72-0012L
.014"	DB72-0014U	DB72-0014L
.016"	DB72-0016U	DB72-0016L
.018"	DB72-0018U	DB72-0018L
.020"	DB72-0020U	DB72-0020L

Thermal

Round	Upper	Lower
.012"	DB73-0012U	DB73-0012L
.014"	DB73-0014U	DB73-0014L
.016"	DB73-0016U	DB73-0016L
.018"	DB73-0018U	DB73-0018L
.020"	DB73-0020U	DB73-0020L

Rectangular	Upper	Lower
.016" x .016"	DB71-1616U	DB71-1616L
.016" x .022"	DB71-1622U	DB71-1622L
.017" x .025"	DB71-1725U	DB71-1725L
.018" x .025"	DB71-1825U	DB71-1825L
.019" x .025"	DB71-1925U	DB71-1925L
.021" x .025"	DB71-2125U	DB71-2125L

Rectangular	Upper	Lower
.016" x .016"	DB72-1616U	DB72-1616L
.016" x .022"	DB72-1622U	DB72-1622L
.017" x .025"	DB72-1725U	DB72-1725L
.018" x .025"	DB72-1825U	DB72-1825L
.019" x .025"	DB72-1925U	DB72-1925L
.021" x .025"	DB72-2125U	DB72-2125L

Rectangular	Upper	Lower
.016" x .016"	DB73-1616U	DB73-1616L
.016" x .022"	DB73-1622U	DB73-1622L
.017" x .025"	DB73-1725U	DB73-1725L
.018" x .025"	DB73-1825U	DB73-1825L
.019" x .025"	DB73-1925U	DB73-1925L
.021" x .025"	DB73-2125U	DB73-2125L

	Stainless Steel	Super Elastic NiTi	Thermal Niti
Round	£24.75	£31.00	£33.75
Rect.	£29.00	£35.00	£38.00

wire & wire products



Super Elastic Nickel Titanium Archwires

EuroLine™ European Arch Shape

This super elastic archwire provides a light and continuous force with a high resistance to deformation. Sold in packs of 10.

Round	Upper	Lower
.012"	DB01-012U	DB01-012L
.014"	DB01-014U	DB01-014L
.016"	DB01-016U	DB01-016L
.018"	DB01-018U	DB01-018L
.020"	DB01-020U	DB01-020L
Rectangular	Upper	Lower
.016" x .016"	DB01-1616U	DB01-1616L
.016" x .022"	DB01-1622U	DB01-1622L
.017" x .025"	DB01-1725U	DB01-1725L
.018" x .025"	DB01-1825U	DB01-1825L
.019" x .025"	DB01-1925U	DB01-1925L
.021" x .025"	DB01-2125U	DB01-2125L
Quantity	Round	Rectangular
1 Pack	£9.50	£11.75
5 Packs	£8.50	£10.75
10 Packs	£8.00	£10.00
20 Packs	£7.50	£9.50

Super Elastic Nickel Titanium on a Reel

Very economical, straightened and wound onto a reel. Sold on 15ft reel.



.012"	DB01-112R	£15.25
.014"	DB01-114R	£15.25
.016"	DB01-116R	£15.25
.018"	DB01-118R	£15.25
Quantity	Price/reel	
5 reels	£14.25	
10 reels	£12.25	

Thermally Active Nickel Titanium Archwires

EuroLine™ European Arch Shape

Thermal wire is an advanced heat activated wire with superior shape/memory properties. This soft wire is easily bent outside the mouth. It gradually increases in force and activates at body temperature ensuring a light, continuous force.

Sold in packs of 10.



Round	Upper	Lower
.012"	DB21-0012U	DB21-0012L
.014"	DB21-0014U	DB21-0014L
.016"	DB21-0016U	DB21-0016L
.018"	DB21-0018U	DB21-0018L
.020"	DB21-0020U	DB21-0020L

Rectangular	Upper	Lower
.016" x .022"	DB21-1622U	DB21-1622L
.017" x .025"	DB21-1725U	DB21-1725L
.018" x .025"	DB21-1825U	DB21-1825L
.019" x .025"	DB21-1925U	DB21-1925L
.021" x .025"	DB21-2125U	DB21-2125L

Quantity	Round	Rectangular
1 Pack	£13.00	£14.85
5 Pack	£11.85	£14.00
10 Pack	£10.60	£12.90
20 Pack	£9.50	£11.75

Beta Titanium

EuroLine™ European Arch Shape

Moderate stiffness for root correction compared to lower modulus materials. Allows for longer activations. Bendability allows for intricate loop designs. Safe for Nickel sensitive patients. Force levels between Nickel Titanium and stainless steel. Sold in packs of 10.

Round	Upper	Lower
.016"	DB23-0016U	DB23-0016L
.018"	DB23-0018U	DB23-0018L

Rectangular	Upper	Lower
.016" x .022"	DB23-1622U	DB23-1622L
.017" x .025"	DB23-1725U	DB23-1725L
.018" x .025"	DB23-1825U	DB23-1825L
.019" x .025"	DB23-1925U	DB23-1925L
.021" x .025"	DB23-2125U	DB23-2125L

Quantity	Round	Rectangular
1 Pack	£29.20	£34.95
5 Pack	£28.50	£31.60
10 Pack	£27.75	£29.75
20 Pack	£27.00	£28.25

wire & wire products

Super Elastic Nickel Titanium Reverse Curve Archwires

Reverse Curve Archwire clinical benefits include:

Ability to open or close patient's bite.
Retraction of flared, protruding incisors.
The prevention of unwanted tooth movement while preserving the arch form.

Sold in packs of 10.



Round	Upper	Lower
.014"	DB51-0014U	DB51-0014L
.016"	DB51-0016U	DB51-0016L
.018"	DB51-0018U	DB51-0018L
Rectangular	Upper	Lower
.016" x .022"	DB51-1622U	DB51-1622L
.017" x .025"	DB51-1725U	DB51-1725L
.018" x .025"	DB51-1825U	DB51-1825L
.019" x .025"	DB51-1925U	DB51-1925L
Quantity	Round	Rectangular
1 Pack	£21.00	£24.50
5 Pack	£20.00	£23.25
10 Pack	£19.00	£21.25
20 Pack	£17.00	£19.10

Stainless Steel Archwires

EuroLine™ European Arch Shape

Gold Tone preformed arches. Round - sold in packs of 100. Rectangular - sold in packs of 10.



Round	Upper	Lower
.012"	DB02-0012U	DB02-0012L
.014"	DB02-0014U	DB02-0014L
.016"	DB02-0016U	DB02-0016L
.018"	DB02-0018U	DB02-0018L
.020"	DB02-0020U	DB02-0020L
Rectangular	Upper	Lower
.016" x .016"	DB02-1616U	DB02-1616L
.016" x .022"	DB02-1622U	DB02-1622L
.017" x .025"	DB02-1725U	DB02-1725L
.018" x .025"	DB02-1825U	DB02-1825L
.019" x .025"	DB02-1925U	DB02-1925L
.021" x .025"	DB02-2125U	DB02-2125L
Quantity	Round	Rectangular
1 Pack	£16.50	£4.25
5 Pack	£15.25	£3.75
10 Pack	£13.50	£3.50
20 Pack	£12.25	£2.75



Stainless Steel Posted Archwires

EuroLine™ European Arch Shape



Gold Tone preformed arches.

Available in 3 sizes: .016 x .022; .018 x .025; .019 x .025

Upper sizes: 34mm, 36mm, 38mm, 40mm

Lower sizes: 24mm, 26mm, 28mm, 30mm

Sold in packs of 10.



Rectangular	Upper	Lower
.016" x .022"	DB02-622UP*	DB02-622LP*
.018" x .025"	DB02-825UP*	DB02-825LP*
.019" x .025"	DB02-925UP*	DB02-925LP*

* Add sizes after UP or LP, e.g. DB02-825UP36 = .018 x .025 36mm

Quantity	Price/pack
1 Pack	£17.50
5 Packs	£16.75
10 Packs	£14.75
20 Packs	£12.50



Archwire Holder

Keeps your work area organised. Stores 12 separate sizes of archwire. Hinged clear lid keeps arches clean. Durable & crystal clear acrylic construction.

DB04-0008

£39.95



Orthogap Disposable Markers

Specially designed thick red wax tip provides an ideal way of marking any wire. Markers sold in packs of 100.


Markers	DB04-0110	£9.50
Dispenser	DB04-0111	£10.25

wire & wire products

Triline™

A new arch shape addressing the deficiencies of existing archwire shape in today's market.

- Natural anterior arch form which helps to ensure longterm stability.
- No over expansion of the Pre-molars.
- Wider in the molar region when compared to standard arch shapes.
- Stainless Steel arches also available with 6 cuspid markings for accurate placement of hooks.
- Sold in packs of 10.



SE Nickel Titanium

	95	100	105
.014	DB31-95/14	DB31-100/14	DB31-105/14
.016	DB31-95/16	DB31-100/16	DB31-105/16
.018	DB31-95/18	DB31-100/18	DB31-105/18
.014 x .025	DB31-95/1425	DB31-100/1425	DB31-105/1425
.016 x .025	DB31-95/1625	DB31-100/1625	DB31-105/1625
.018 x .025	DB31-95/1825	DB31-100/1825	DB31-105/1825
.019 x .025	DB31-95/1925	DB31-100/1925	DB31-105/1925
.021 x .025	DB31-95/2125	DB31-100/2125	DB31-105/2125

Thermal Nickel Titanium

	95	100	105
.014	DB33-95/14	DB33-100/14	DB33-105/14
.016	DB33-95/16	DB33-100/16	DB33-105/16
.018	DB33-95/18	DB33-100/18	DB33-105/18
.014 x .025	DB33-95/1425	DB33-100/1425	DB33-105/1425
.016 x .025	DB33-95/1625	DB33-100/1625	DB33-105/1625
.018 x .025	DB33-95/1825	DB33-100/1825	DB33-105/1825
.019 x .025	DB33-95/1925	DB33-100/1925	DB33-105/1925
.021 x .025	DB33-95/2125	DB33-100/2125	DB33-105/2125

Stainless Steel

	95	100	105
.019 x .025	DB32-95/1925	DB32-100/1925	DB32-105/1925
.021 x .025	DB32-95/2125	DB32-100/2125	DB32-105/2125

Stainless Steel with Cuspid Markings

	95	100	105
.019 x .025	DB32-95/1925M	DB32-100/1925M	DB32-105/1925M
.021 x .025	DB32-95/2125M	DB32-100/2125M	DB32-105/2125M

	Round	Rectangular
SE Nickel Titanium	£12.00	£16.00
Thermal Nickel Titanium	£17.00	£21.50
Stainless Steel	£7.50	
Stainless Steel with Cuspid Markings	£21.50	



wire & wire products



A.J. Wilcock® Australian Wire

Although Australian wire has been imitated by many manufacturers over the years, A.J. Wilcock Scientific & Engineering remains the global market leader in high tensile stainless steel archwires.

For over 50 years, "Australian Wire" has been a household name in orthodontics.

The extraordinary properties of these wires are due to the research and fine engineering skills employed in the processing of the wires. Wilcock wires are well known for their resilience and ability to withstand masticatory forces as well as being able to maintain their shape even when auxiliaries and elastics are used. There is no other wire which opens the bite as effectively as Wilcock wire.

The wires on this page are listed in order of resilience.

The Regular and Regular Plus are easily formed and are excellent wires for general use and utility wire. The Special and Special Plus wires are ideal for bite opening and where high resilience is required.

Australian Wire - Regular Grade

	25ft Coil	Price	10" Lengths*	Price
.012"	DB02-229-000	£13.00		
.014"	DB02-230-100	£13.00	DB02-230-101	£19.00
.016"	DB02-231-300	£13.00	DB02-231-301	£19.00
.018"	DB02-232-400	£13.00	DB02-232-401	£19.00
.020"	DB02-233-500	£13.00	DB02-233-501	£19.00

Australian Wire - Regular Plus Grade

	25ft Coil	Price	10" Lengths*	Price
.014"	DB02-230-110	£13.00	DB02-230-111	£20.00
.016"	DB02-231-310	£13.00	DB02-231-311	£20.00
.018"	DB02-232-410	£13.00	DB02-232-411	£20.00
.020"	DB02-233-510	£13.00	DB02-233-511	£20.00

Australian Wire - Special Grade

	25ft Coil	Price	10" Lengths*	Price
.014"	DB02-230-120	£13.00	DB02-230-121	£22.75
.016"	DB02-231-320	£13.00	DB02-231-321	£22.75
.018"	DB02-232-420	£13.00	DB02-232-421	£22.75
.020"	DB02-233-520	£13.00	DB02-233-521	£22.75


Australian Wire - Special Plus Grade

	25ft Coil	Price	10" Lengths*	Price
.014"	DB02-230-130	£19.00	DB02-230-131	£29.00
.016"	DB02-231-330	£19.00	DB02-231-331	£29.00
.018"	DB02-232-430	£19.00	DB02-232-431	£29.00
.020"	DB02-233-530	£19.00	DB02-233-531	£29.00

*The 10" lengths are supplied in packs of 10.

Lingual ArchWires

DB Orthodontics Lingual Archwires have been developed to meet the challenging mechanics of modern lingual treatment. Archforms are designed and adapted precisely to lingual anatomy. Various alloys are offered to maximise treatment efficiencies in the narrow interbracket dimensions of lingual orthodontic treatment.



A.J. Wilcock® Australian Wire Stainless Steel Lingual Arches


A.J. Wilcock Combination Archwires are rectangular in the anterior segment to maintain torque while the round posterior segments allow sliding mechanics with reduced friction.

Size	Upper	Price
Anterior .016" x .022" Posterior .016"	DB14-1622/16	£29.00/pk 2
Anterior .017" x .025" Posterior .017"	DB14-1725/17	£29.00/pk 2
Anterior .019" x .025" Posterior .019"	DB14-1925/19	£29.00/pk 2
Anterior .021" x .025" Posterior .021"	DB14-2125/21	£29.00/pk 2
Anterior .021" x .027" Posterior .021"	DB14-2127/21	£29.00/pk 2

Stainless Steel

available in .014", .016", .018", .016" x .022"

Size	Upper	Price
.014"	DB14-0014	£13.00/pk 5
.016"	DB14-0016	£13.00/pk 5
.018"	DB14-0018	£13.00/pk 5
.016" x .022"	DB14-1622	£13.00/pk 5



Beta Titanium


available in .016", .016" x .022", .0175" x .0175"

Size	Upper	Price
.016"	DB14-0016B	£37.00/pk 5
.016" x .022"	DB14-1622B	£37.00/pk 5
.017" x .025"	DB14-1725B	£37.00/pk 5

wire & wire products

Lingual NiTi Archforms

(arches shown actual size)



Span of anterior segment

Size 1 Upper	32mm	Size 1 Lower	27mm
Size 2 Upper	33.5mm	Size 2 Lower	30mm
Size 3 Upper	36mm	Size 3 Lower	31mm

Nickel Titanium

available in .010", .012", .014", .016"

Size		Upper	Lower	Price
Size 2	.010"	DB13-S2/010U	DB13-S2/010L	£9.25/pk 2
Size 2	.013"	DB13-S2/013U	DB13-S2/013L	£9.25/pk 2
Size 2	.016"	DB13-S2/016U	DB13-S2/016L	£9.25/pk 2


* Add sizes after S, e.g. DB13-S1/010U = .010 32mm

Thermal Nickel Titanium

available in .016" .016" x .022", .017" x .017", .017" x .025"

Size		Upper	Lower	Price
Size 2	.016"	DB13-S2/016TU	DB13-S2/016TL	£9.75/pk 2
Size 2	.016" x .022"	DB13-S2/1622U	DB13-S2/1622U	£9.75/pk 2
Size 2	.017" x .017"	DB13-S2/1717U	DB13-S2/1717U	£9.75/pk 2
Size 2	.017" x .025"	DB13-S2/1725U	DB13-S2/1725U	£9.75/pk 2

* Add sizes after S, e.g. DB13-S1/1622U = .016" x .022" 32mm



Braided Archwires

Lower force than solid stainless steel archwires. Sold in packs of 10.

.016" x .022"	DB62-1622	£15.25
.017" x .025"	DB62-1725	£15.25
.019" x .025"	DB62-1925	£15.25
.021" x .025"	DB62-2125	£15.25



Co-Ax Stainless Steel on a Spool

This Co-Ax design of five wires wrapped tightly around a single core wire produces a resilient wire that is suitable for the initial arch stage of treatment. Sold on 30ft spool.

.0155	DB62-0155	£14.25
.0175	DB62-0175	£14.25



NiTi Open Coil Spring on a Spool

For opening and closing spaces with a gentle and constant force. Cut to length required, fits over any archwire. Sold on 15" spool.

Light .009" x .030"	DB41-0930	£13.75
Medium .010" x .036"	DB41-1036	£13.75
Heavy .012" x .036"	DB41-1236	£13.75



NiTi Closed Coil Springs on a Spool

Cut to length required, fits over any archwire. .10" x .036" Medium. Sold on 15" Spool.

Closed	DB41-1036C	£16.00
--------	------------	---------------



Stainless Steel Coil Springs

Stainless steel coil springs are precision wound in two styles - open and closed. 3ft spool. Wire diameter .010 (.25mm) x .030 (.76mm)

Closed	DB02-1030C	£7.50
Open	DB02-1030	£7.50



NiTi Preformed Springs

Nickel Titanium Preformed Springs close or open spaces at a level which is unapproachable with stainless steel springs. Sold in packs of 10.

Force	9mm	12mm	Price/pack
Light 150g	DB41-0150	DB41-0151	£27.00
Medium 200g	DB41-0200	DB41-0201	£27.00
Heavy 250g	DB41-0250	DB41-0251	£27.00

wire & wire products



Preformed Coated Ligatures

.012" Pk 100
DB02-0212LIG
£15.00




Pre-cut Kobayashi Short Ligatures

.012" Pk 100
DB02-112LIG
£23.00



Pre-cut Short Ligatures

Pk 1000
.008" DB02-108LIG
.010" DB02-110LIG
£21.00



Pre-cut Short Coated Ligatures

.012" Pk 100
DB02-0312LIG
£18.25



Preformed Ligatures


Pk 1000
.008" DB02-008LIG
.010" DB02-010LIG
£ 13.50



Preformed Kobayashi Ligatures

Pk 100

.012"	DB02-012LIG	£7.50
.014"	DB02-014LIG	£7.50



NiTi Ice Instrument

Turns NiTi wire soft and flexible immediately.

DB05-0322
£16.00



NiTi Molar Distalizing Spring

Distalize molars without headgear or lip bumper. Fits over any size archwire. .010 x .045 open/wound. Sold in packs of 3 x 7" lengths.

DB41-1045
£27.25



Mathieu Needle Holder

DB05-0406
£37.50



Straight Smaha

Very fine tip.
DB05-0409
£17.00