

# create™

orthodontic alginate  
& accessories



Natural fruit aroma

Attractive Fuchsia coloured powder

Dust Free

Extra fast setting

*See over for more details.*

freephone \* 0800 783 3552

# create™

## orthodontic alginate



**Create™** Orthodontic Alginate has been specifically formulated for Orthodontic use. Create™ has an extra fast setting time, high elasticity and is tear resistant.

**Create™** Orthodontic Alginate is an attractive fuchsia coloured powder with a red fruit flavour.

**Create™** is also dust free: The powder is treated in such a way that no dust is generated during dosing and mixing. This facilitates the mixing process through easy absorption of the water by the powder, offering the orthodontist convenient and hygienic handling.

**Create™** Orthodontic Alginate is sold in packs of 500g.

	No of packs	Total cost	Cost/500g
500g	1	<b>£8.85</b>	<b>£7.60</b>
5kg	10	<b>£72.50</b>	<b>£7.25</b>
10kg	20	<b>£135.00</b>	<b>£6.75</b>
20kg	40	<b>£254.00</b>	<b>£6.35</b>



**Alginate Container**

DB04-0276  
**£31.30**



**Water Dosing Bottle**

DB04-0275  
**£18.30**



6

**freephone \* 0800 783 3552**

# Tropical Mixing Bowls

Beautiful transparent colours. Highly flexible to ensure smooth mixing.

Colour	Bowls	Price	Spatulas	Price
Yellow	DBL5-0003	£6.50	DBL5-0032	£4.00
Pink	DBL5-0004	£6.50	DBL5-0033	£4.00
Green	DBL5-0005	£6.50	DBL5-0034	£4.00
Orange	DBL5-0006	£6.50	DBL5-0035	£4.00
Blue (not shown)	DBL5-0007	£6.50	DBL5-0036	£4.00



# Lock Tight Impression Trays

Designed specifically for Orthodontists. Retention slots and extra perforations provide superior retention. Colour & number coded for easy identification.

Colour	Upper	Lower
1 Red	DB04-0230U	DB04-0230L
2 Orange	DB04-0231U	DB04-0231L
3 Blue	DB04-0232U	DB04-0232L
4 Green	DB04-0233U	DB04-0233L
5 Yellow	DB04-0234U	DB04-0234L
6 Purple	DB04-0235U	DB04-0235L

Quantity Rates	Price/pack
1 pack (50 trays)	£18.50
5 packs	£16.50
10 packs or more	£14.25



## Imp-Lock Trays

Durable plastic, excellent alginate retention.

Size	Colour	Upper	Lower
Xtra Small	Red	DB04-0215U	DB04-0215L
Small	Orange	DB04-0214U	DB04-0214L
Medium	Blue	DB04-0213U	DB04-0213L
Large	Green	DB04-0212U	DB04-0212L
Xtra Large	Yellow	DB04-0211U	DB04-0211L

Quantity Rates	Price/pack
1 tray	£0.70
50 trays	£0.65
100 trays	£0.55
300 trays	£0.50
<b>Starter Kit (50)</b> 5 of each size DB04-0210S/Kit	<b>£29.00 = 0.56/each</b>

**db** \*  
orthodontics

**DB Orthodontics Limited**  
Ryefield Way Silsden West Yorkshire  
BD20 0EF United Kingdom

**telephone** +44[0] 1535 656 999  
**fax** +44[0] 1535 656 969  
**email** sales@dbortho.com  
**internet** www.dbortho.com

## Rensa™

### Impression Disinfectant Powder 700 grams

Fast disinfectant for all kinds of impression material such as silicone, alginates and polyether.

Bactericidal in 5 minutes.

Contents:

Potassium peroxomonosulphate 50%,  
Sulphamic acid 5%, Sodium alkyl benzene 15%.

4509-0070

**£18.50**



### Impression Disinfectant Bath

3 litre instrument bath. Inner removable draining tank with lid.

4509-0071

**£28.95**



## Model Boxes

		Price
8" x 3" x 3"	DB04-0550	£2.00
9" x 3" x 3"	DB04-0551	£2.00
11" x 3" x 3"	DB04-0552	£2.00

freephone \* **0800 783 3552**