

interdental stripping

Orthospace - Interproximal Reduction System

for safe & fast interproximal reduction

New and improved interproximal stripping treatment

The Orthospace system is a thin stainless steel flexible interproximal file, which facilitates the contouring of enamel and will reduce many types of dental materials. It enables the often necessary reduction of the interproximal tooth structure, correcting inadequate space, caused by an insufficient arch length. Interproximal re-contouring is easily achieved in a minimum time with a more natural finish.

A safer, faster, more precise and reliable solution

Compared to the other manual systems or with discs, the Orthospace files combined with the handpiece (EVA system), are a more efficient, faster and practical instrument. In fact this system allows more freedom of movement during the treatment. Thanks to the oscillating movement, a safe treatment is guaranteed to the patient.

Using Orthospace removal of excess tooth structure up to 14mm per dental arch is achievable and borderline extraction cases are avoided. Single sided files are perfect for situations where you want to reduce tooth structure on one side only without any damage or marking to the adjacent tooth. Suitable for cases where you may have a crown on an adjacent tooth which you want to keep intact.

The Kit includes:

- 1x Eva Handpiece
- 1x Extractor
- 1x Pack of files 7pcs
- 1x Set of wedges
- 1x Tweezers for wedges
- 1x Set of measuring gauges

Kit DB06-DC-IRS **£61.00**

Replacement Strips Pk 3 **£56.50**



The Orthospace System can be used during treatment.

	Left Sided Single Pack	Right Sided Single Pack
Grain	Product Code	Product Code
90	DB06-OC-L90-B1	DB06-OC-R90-B1
60	DB06-OC-L60-B1	DB06-OC-R60-B1
40	DB06-OC-L40-B1	DB06-OC-R40-B1
25	DB06-OC-L25-B1	DB06-OC-R25-B1
15	DB06-OC-L15-B1	DB06-OC-R15-B1
	Left Sided Pack 3	Right Sided Pack 3
90	DB06-OC-L90-B3	DB06-OC-R90-B3
60	DB06-OC-L60-B3	DB06-OC-R60-B3
40	DB06-OC-L40-B3	DB06-OC-R40-B3
25	DB06-OC-L25-B3	DB06-OC-R25-B3
15	DB06-OC-L15-B3	DB06-OC-R15-B3


Double sided strips are available in all grains, as well as 150 and 120 grains.



Stripping Discs

Serrated diamond discs with fine diamond grit, for separating ceramics, plaster, acrylics, resin veneers. Available in 2 diameter sizes: 19mm (20,000 rpm) and 22mm (15,000 rpm).

DB04-1000	19mm Double sided	£14.25/each
DB04-1001	22mm Double sided	£14.25/each



Superflex Discs

Double sided diamond discs with fine diamond grit, for separating ceramics, plaster, acrylics and resin veneers. Available in three diameter sizes. Fixed mandrel.

DB04-1010	16mm Double sided	£10.25/each
DB04-1011	19mm Double sided	£10.25/each
DB04-1012	22mm Double sided	£10.25/each
DB04-1013	19mm Single sided Upper	£10.25/each
DB04-1014	22mm Single sided Upper	£10.25/each
DB04-1015	19mm Single sided Lower	£10.25/each
DB04-1016	22mm Single sided Lower	£10.25/each



Disc Guards

For safe, efficient and effective reproximation. Protects patient's lips, cheeks and gums from laceration. Easy to use: attach discguard to adaptor, slide adaptor onto slow speed straight handpiece, insert disc and close chuck latch.

DB04-1002
£18.50/set of 3



Interproximal Diamond Bur

For interproximal stripping.

DB04-1017	Medium - 524	£19.00/pack 5
DB04-1018	Super Fine - 504	£19.00/pack 5
DB04-1019	Fine - 514	£19.00/pack 5



Adaptor

Fits most slow speed straight handpieces.

DB04-1003
£20.75



Interproximal Stripper Kit

For removal & recontouring of interproximal enamel. Kit comes with six 4mm double sided abrasive strips, medium grit. Replacement strips sold in packs of 12.

DB05-0324	Enamel stripper kit	£28.75
DB05-0324B	Abrasive strips double sided 4mm	£11.25/pack 12
DB05-0324S	Abrasive strips single sided 4mm	£9.75/pack 12