


# headgear & accessories



## High Pull Headgear

Maximum comfort washable cotton. Sold as each.

Green	●	DB08-0050	£5.75
Denim	●	DB08-0051	£5.75
Pink	●	DB08-0052	£5.75
Blue	●	DB08-0053	£5.75
Purple	●	DB08-0054	£5.75
Black	●	DB08-0055	£5.75



## Neck Pads

Maximum comfort washable cotton. Sold as each.


Green	●	DB08-0034	£4.50
Denim	●	DB08-0035	£4.50
Pink	●	DB08-0036	£4.50
Blue	●	DB08-0037	£4.50
Purple	●	DB08-0038	£4.50
Black	●	DB08-0039	£4.50



## Masel Safety Straps

Convert regular facebows into safety Facebows with these inexpensive straps. Sold in packs of 10.

DB80-0165  
£13.50



## Safety Release Modules

Safety of the system is provided by the clips which release the spring when traction force is higher than 2 kilos each side. This system reduces the hazard of slingshot rebound of facebow in case of disconnection from headgear tubes. Coil springs are incorporated into the modules to provide constant force. Available in 3 force levels. Sold in packs of 10.

Light Clear 600g	LM0810-01	£30.50
Medium Yellow 750g	LM0810-02	£30.50
Heavy Blue 1000g	LM0810-03	£30.50




### Extra Oral Facebows

Long Outer bow 1.80mm with looped Inner bow 1.14mm (.045") in diameter. Laser welded joint with bright electro polishing for maximum hygiene.

*Caution: All Extra Oral Facebows must be used with safety release modules only.*

Long Outer Bow 83mm	LM4805-01	<b>£7.30</b>
Long Outer Bow 90mm	LM4805-02	<b>£7.30</b>
Long Outer Bow 97mm	LM4805-03	<b>£7.30</b>
Long Outer Bow 104mm	LM4805-04	<b>£7.30</b>
Long Outer Bow 111mm	LM4805-05	<b>£7.30</b>



### Short Outer Facebows

Inner bow .045". Short Outer bow.

Short Outer Bow 83mm	LM4835-01	<b>£7.25</b>
Short Outer Bow 90mm	LM4835-02	<b>£7.25</b>
Short Outer Bow 97mm	LM4835-03	<b>£7.25</b>
Short Outer Bow 104mm	LM4835-04	<b>£7.25</b>
Short Outer Bow 111mm	LM4835-05	<b>£7.25</b>



### Multi Adjustable Face Mask

Fully adjustable to any patient; lightweight design. Two replacement pads included.

Blue Face Mask	DB08-0060	<b>£38.00</b>
Red Face Mask	DB08-0061	<b>£38.00</b>



### Headgear Bags

Sold in packs of 10.

LM3620-01  
**£19.50**