

instruments



Mosquito Forceps

High quality instrument. Fine serrated tips for placing elastomeric ligatures. Matt finish. Autoclavable.

DB05-0315

£8.50



Mosquito Forceps with Grooved Tip

Grooved tip means a non-slip grip when placing elastic ligatures.

DB05-5316

£13.25



Mosquito Forceps with Curved Tip

Same high quality as standard forceps. Curved tip. New matt finish. Autoclavable.

DB05-0315C

£11.50




Mosquito Forceps with Hook

Same quality as standard mosquito. Hook to stop ligatures slipping.

DB05-0315H

£9.75



Direct Bond Holders

Style A Deluxe Holder

Fine tip with positioning blade. Reverse action holds all types of brackets and releases smoothly.

DB05-0421 £11.50

Style B Standard Holder

Standard holder with positioning blade. Serrated tips. Reverse action for positive holding.

DB05-0420 £10.25

Style C Curved Holder

DB05-0418 £11.50



Reverse Keat Bondable Tube Holder

Helps with accurate positioning of bondable tubes.

DB05-0423

£22.50



Positive Positioner

For precise placement of posterior direct bond attachments. Lightweight design, spring release action.

DB05-0323

£26.00



College Tweezers

High quality with serrated tips. New matt finish.

DB05-0319

£7.75



NiTi Ice Instrument

Turns NiTi wire soft and flexible immediately. Fill the NiTi Ice Instrument with water and freeze. The lightweight compact design allows for pinpoint activation of nickel titanium archwires. Wire turns soft and flexible immediately. Works well in anterior or posterior regions.

DB05-0322

£16.00



Surgery Adams Plier

Tungsten Carbide Tips. Ideal for all orthodontic appliance construction and adjustment.

DB05-0430

£41.00



Large Adams Plier

Tungsten Carbide Tips. Ideal for all orthodontic appliance construction and adjustment.

DB05-0430L

£45.00

instruments



Adams No.64

Square 1mm tips are parallel at 1mm opening. Fine Tungsten Carbide coated impregnated tips.

DB05-0064

£90.00



Spring Former No.65

Fine Tungsten Carbide Impregnated coating to ensure excellent grip on all sizes of wire.

DB05-0065

£90.00



Pedo Curved Scissors

DB05-0331

£12.50



Beebee Straight Scissors 10cm

DB05-0332

£12.50



Crown & Gold Scissors 12cm

DB05-0333

£12.50



Boon Gauge

Helps with precise bracket placement. Accurately measures height of bracket from edge. 3.5, 4, 4.5, 5mm. Made from surgical stainless steel. Autoclavable.

DB05-0004

£12.50



Ruler

150mm 6" Stainless Steel ruler. Autoclavable.

DB05-0318

£8.50



Open Bite/Wedge Ruler

Ideal for diastemes, expansion appliances. Measuring range 0.4mm - 8mm. Autoclavable.

DB05-0329

£9.25



Divider

Overall length 4" (115mm). Measuring range up to 95mm. Can be used to measure models and intraorally. Made from surgical stainless steel. Autoclavable.

DB05-0

£13.75




Force Module Separator

Easy spring-back action. Undercut tips positively hold separation modules.

DB05-0302

£61.75

instruments



Autoclavable Heavy Duty Cutter with TC Tips

DB05-0500T
£36.75



Autoclavable Bite Sticks

Tungsten Carbide Tips. Ideal for all orthodontic appliance construction and adjustment.

Red	DB05-5100	£6.25
Light Green	DB05-5101	£6.25
Purple	DB05-5102	£6.25
Turquoise	DB05-5103	£6.25
Orange	DB05-5104	£6.25
Cobalt Blue	DB05-5105	£6.25



High Heat Bite Sticks

A major improvement to the traditional orthodontic bite stick design. High Heat Bite Stick is manufactured from a high impact plastic that is similar to stainless steel and will not flex or bend like older bite sticks. A sure grip handle is provided that makes it much easier to grip, even when the surface is wet. Withstands up to 600°F. Guaranteed not to melt.

DB05-5053
£20.00



Mirrors No.4

Tested to withstand 2000 cycles in an autoclave. The design of the reflecting layers provides maximum image quality and extremely high resistance to deterioration. Heads packed in boxes of 12. Handles packed individually.

Mirror Heads	DB05-0313M	£13.50/Pk12
Handles	DB05-0313	£2.00
Mirror Handle & Head	DB05-0313H	£3.00



Mathieu Needle Holder

High quality holder with T.C. tips for placement of stainless steel ligatures.

DB05-0406
£37.50



Grooved Needle Holder

Grooved tip.
DB05-0408
£17.75



Straight Smaha

Very fine tip.
DB05-0409
£17.00



Short Probe

DB05-0304
£5.75



Distal Bender

Works on all stainless steel, most NiTi and Beta III wires.
Can be used on all sizes of wires up to .021" x .028".
DB05-0321
£17.50



Double Ended Tucker

The most popular ligature director. One end notched, one end curved. For directing archwires and ligatures.
DB05-0311
£12.50

instruments



Cement Spatula

DB05-0325
£10.00



Twirl-On

DB05-5317
£11.25



Probe

DB05-0312
£11.75



Flat Plastic No.6

DB05-0316
£8.50



Mitchell Trimmer No.4

DB05-0317
£11.00



Wards Carver No.2

DB05-0320
£9.50

DBO instruments

- > Diamond honed edges.
- > Precision tip alignment.
- > Withstand all accepted methods of sterilisation.
- > Manufactured from Surgical Grade Stainless Steel.
- > Come with a 5 year warranty against corrosion.



Small Distal End Cutter

Small head provides access to hard to reach areas while minimising cheek impingement.

Maximum cutting capacity: .018" x .025"

Minimum cutting capacity: .014"

DB05-0205S

£56.50



Long Handled Distal End Cutter

Small head on this cutter provides safety hold features and long handles for that extra reach.

Maximum cutting capacity: .018" x .025"

Minimum cutting capacity: .014"

DB05-0205L

£56.50



Hard Wire Cutter

Standard wire cutter, non-angled, designed to cut hard wires. Compact design. Stainless steel with hard chrome plating.

Maximum cutting capacity: .022" x .028".

DB05-0209

£56.50



Hard Wire Cutter 15°

Specially designed for cutting hard wires. The cutting angle is 15° and the cutting edge has been diamond honed.

Maximum cutting capacity: .022" x .028".

DB05-0215

£56.50



Black Hard Wire Cutter 15°

Distinctive black oxide finish. Features extra strong cutting edges and a 15° cutting angle.

Cuts archwires up to .022" x .028".

Autoclaving is not recommended.

DB05-0215B

£56.50

instruments



Mini Pin & Ligature Cutter

This cutting plier has finer tips for easy access into difficult areas. Cuts pins, ligatures and elastics. Soft wire up to .015"

DB05-0203

£56.50



Weingart Plier

The most popular utility plier. Serrated tips hold wire firmly at a convenient working angle. Tapered tip fits easily between brackets. Rounded for safety.

DB05-0100

£56.50



Bird Beak Plier

Popular and versatile utility plier for working round wire up to .030" in diameter. Tips are designed to prevent scoring.

DB05-0101

£56.50



Three Jaw Plier

Precision-aligned tips for consistent bends. Gradually tapered for fine wire bending.

Round edges for wire up to .030".

DB05-0103

£56.50



Posterior Band Remover

The nylon tip rests on the occlusal of the molar and the metal tip catches the edge of the band for upward or downward removal. Available long and short.

Posterior Band Remover Long	DB05-0102	£56.50
Posterior Band Remover Short	DB05-0115	£56.50



Direct Bond Bracket Remover

The bond remover tips wedge between both edges of the base of the bracket and the tooth surface. Lifting off the bracket.

DB05-0112

£56.50



Angled Bracket Remover

Designed to remove all bonded attachments. The 60° angled tip allows for easy access to difficult areas.

DB05-0344
£56.50



Bracket Remover with Pad

Removes all direct bond brackets and attachments. Nylon tip rests on incisal/occlusal edge of teeth while the metal edge fits under bracket flange.

Bracket Remover with Pad	DB05-0113	£56.50
Replacement Delrin Pads	DB05-0308	£20.50



Crimpable Archwire Hook Plier

Designed for surgical ball hooks and auxiliaries, the slotted tip holds auxiliaries secure to allow for tight crimping.

DB05-0172
£56.50



Arch Former - Tweed Style

Tweed style plier for square or rectangular wire. Edges will not score wire. Will bend wires up to .022" x .028".

DB05-0104
£56.50



Spring Former Plier

Designed for round wire up to .016" in diameter. The longer tapered beaks make it easier to bend small diameter loops. Prevents wire scoring.

DB05-0106
£56.50



Nance Plier

Precision loop forming with 4 steps tip 3, 4, 5 and 6mm. Excellent for wires up to .022". Bevelled edges prevent wire scoring and one serrated beak prevents wire slippage.

DB05-0116
£56.50

instruments



Omega Loop Former

Concave beak is contoured for making more uniform Omega loops. 3 sections of .045", .060" and .075" in diameter and .090" in length.

DB05-0108

£56.50



Wire Contouring Plier

The smooth inner surface allows contouring and forming of archwires. Forms without scoring or torquing. Perfect for accentuating curves.

DB05-0110

£56.50



Jarabak Plier

Precision grooves assure accurate bending and closing of loops. For wires up to .020".

DB05-0107

£56.50




Straight Utility - How Style

A general utility plier with 3/32" diameter tips which are serrated. Designed with a rounded tip for maximum patient comfort and safety in mind.

DB05-0105

£56.50




Utility Plier Angled 15°

A general utility plier, the same as 0105 with machined tips, 3/32" in diameter, but at a 45° offset angle. Serrated tips hold wires firmly.

DB05-0109

£56.50



Turrets for Forming Rectangular or Round Archwires

For forming archwires with or without torque adjustment. Available in 3 different versions & colours. The wire slotted body is made of stainless steel.

Blue Turret

For forming rectangular archwires .016" up to .0215".


Gold Turret

For forming rectangular archwires with torque adjustment: 0° -7° -10° -13° -16°. For use only with wire .018".

Silver Turret

For forming rectangular archwires with torque adjustment: 0° -7° -10° -13° -16°. For use only with wire .022".

Blue Turret	LP1167-00	£101.00
Gold Turret	LP1168-00	£101.00
Silver Turret	LP1170-00	£101.00



Instrument Code Rings

Designed to fit orthodontic plier handles, scalers, band seaters, hemostats and other orthodontic instruments.

Made out of high temperature silicone that can withstand repeated autoclaving and heat sterilisations up to 370 degrees.

Available in packs of 100.

£8.75

Blue	ix110	Orange	ix115	Yellow	ix119
Green	ix111	Pink	ix116	Brown	ix120
Purple	ix112	White	ix117	Dark Grey	ix121
Red	ix113	Tan	ix118	Mixed Colours	ix122
Black	ix114				



Replacement Tips

Pack comes with 2 tips, screws and a wrench.

Metal Tips Wide	DB05-0309	£20.50
Metal Tips Narrow	DB05-0310	£20.50



Delrin Pads

Heat sterilisable replacement pads. Sold in pack of 4.

1/4" pads	DB05-0307	£20.50
3/16" pads	DB05-0308	£20.50