


photographic supplies



Rhodium Photographic Mirrors

High resolution intra-oral photos depend on superior reflective mirrors. We offer four styles that give you unparalleled access to all intra-oral regions. Tough chromium mirror surface is permanently bonded to a metal base that resists chipping. Rhodium intra-oral mirrors are the best choice for photography. Coated both sides for durability. Rhodium provides higher reflectance for sharper, precise images with 35mm and digital photography. Rhodium mirrors withstand repeated steam and autoclave sterilization. Sterilizable up to 350° F, 180°C. Sold as each.

A Palatal Adult	DB04-0160	£40.15
B Palatal Child	DB04-0161	£40.15
C Buccal	DB04-0162	£40.15
D Lingual	DB04-0163	£40.15



Rhodium Mirror Handle

Mirror Handle for use with rhodium mirrors, make taking photos much easier. Fully autoclavable.


Rhodium Mirror Handle	DB05-0326	£70.00
Rhodium Mirror Handle with Mirror		£99.00



Metallic Blades for Photo

2 sizes: small (no.1) & large (no.2). Mirror polished effect. For intra-oral photography. Their specific shape allows very good results. For occlusal photos use large. For lingual photos use small. Sold as each.

Small	DB04-0165/1	£40.55
Large	DB04-0165/2	£40.55



Rhodium Mirrors with Handle

Same high quality as standard mirrors, with the addition of a S/S handle. Available in Buccal (1) and Occlusal (2). Sold as each.

Buccal (1)	DB04-0159	£88.05
Occlusal (2)	DB04-0158	£88.05



Parallelism Mirror

High quality mirror with serigraphy allowing to check parallelism. Sold as each.

Parallelism Mirror	DB04-0169	£34.80
--------------------	-----------	--------



Large Mirrors

Special for intra-oral photography. Complete with handle. Sold as each.

Mirror ø 4cm	DB04-0166	£26.55
Mirror ø 5cm	DB04-0167	£34.80



Cheek Retractors (autoclavable)

Very clear, excellent for photography. Can be autoclaved at 134°C or cold sterilised. Sold in pairs.

Small	DB04-0174	£26.00
Large	DB04-0175	£26.00



Double Ended Photo Retractors (autoclavable)

Sized for adult and child. Durable plastic design. Sold in pairs.

Double Ended Photo Retractors	DB04-0183	£14.95
-------------------------------	-----------	---------------



Photo Retractors (autoclavable)

A photo retractor design that provides a buccal extension to hold the lips and cheeks away from the teeth, resulting in improved intraoral photography. Sold in pairs.


Small	DB04-0176	£10.95
Large	DB04-0177	£10.95



Photo Cheek Retractors (autoclavable)

Each end of the retractor has a unique size that can be used with children and adults. The ends are elongated, unlike many other photo retractors resulting in superior intra oral photography. The photographer is able to take pictures without the aid of an assistant. Sold in pairs.


Small	DB04-0178	£25.60
Large	DB04-0178L	£25.60



Kombee (autoclavable)

Double ended occlusal and buccal cheek retractor. Keeps lips retracted for better occlusal pictures. Sold as each.

Small	DB04-0189	£9.00
Large	DB04-0188	£9.00



Cheek Retractors (autoclavable)

Long handled. Sold in pairs.

Small	DB04-0190	£10.00
Large	DB04-0191	£10.00

photographic supplies



Hand-held Cheek Retractors (autoclavable)

Smooth, rounded edges. Hand-held for increased visibility/accessibility. Fully autoclavable. One universal size fits all patients. Ideal for intra-oral photography.

Cheek Retractors	DB04-0184	£11.85
------------------	-----------	---------------



Soft Spandex (autoclavable)

Soft material allows greater flexibility. Soft cheek retractors are less prone to fracture than hard plastic type. Ideal for intra-oral photography as matt finish does not reflect glare from flash unit. Fully autoclavable. Standard and junior size.

Regular	DB04-0185	£17.30
Small	DB04-0186	£17.30



Extra Oral Cheek Expander (autoclavable)

Autoclavable Extra-oral Cheek Expanders are specifically designed to resist fracture. The unique tab allows the cheeks to be pulled out further. Made from a heavy duty plastic, it is autoclavable without risk of cracking or breaking. Sold in pairs.


Regular	DB04-0170	£12.45
Small	DB04-0170S	£12.45



Lip Ring Retractor

The spring like action maintains the retraction of the patient's lips. Sized for adult and children. Ideal for bonding and intra-oral photography. Cold sterilizable.

Lip Ring Retractor	DB04-0187	£9.90
--------------------	-----------	--------------



Retract Ezz Photo Lip Retractor

Innovative design provides maximum exposure to occlusal photographs. Can be cold sterilised or heat sterilised up to 290°F/143°C. Sold as each.


Retract Ezz Retractor	DB04-0154	£19.50
-----------------------	-----------	---------------



Contrastors (autoclavable)

Black anodized aluminium. Suitable for normal and digital photography.
 Cleaning: only water with mild soap. Sterilisation: autoclave to 134°C. Gives contrast to your images and eliminates useless elements. Sold as a set of 2 pieces: 1 wide + 1 narrow.

Contrastors	DB04-0157	£38.05
-------------	-----------	---------------



Sternberg Lip Retractors (autoclavable)

Sold in pairs.

A		
120 mm (4 6/8 ")	DB04-0193	£11.20
140 mm (5 1/2 ")	DB04-0194	£11.90
160 mm (6 2/8 ")	DB04-0195	£12.45
B		
140 mm (5 1/2 ")	DB04-0196	£11.90
160 mm (6 2/8 ")	DB04-0197	£12.45
180 mm (7 1/16 ")	DB04-0198	£13.10



Stainless Steel Cheek Retractor (autoclavable)

Glare free brushed metal excellent for photography.
 Sold as each.

Cheek Retractor	DB04-0199	£19.95
-----------------	-----------	---------------

Photographic Sterilisation Cassettes

See page 101 for full details.

